

# Certificate of Compliance

Certificate Number 090905 - E215059  
Report Reference E215059, June 17th, 2005  
Issue Date 2005 September 9

Page 1 of 2



*Issued to:* **EPPENDORF A G**  
**BARKHAUSENWEG 1**  
**22339 HAMBURG GERMANY**


*This is to certify that representative samples of* **Centrifuge**  
**Model 5418**

*Have been investigated by Underwriters Laboratories Inc.® in accordance with the Standard(s) indicated on this Certificate.*

*Standard(s) for Safety:* See Addendum for Standards

*Additional Information:* ELECTRICAL RATING:  
Voltage: 120 V ac  
Frequency: 50-60 Hz  
Current: 2.6 A  
Power: 170 W

**Only those products bearing the UL Listing Mark for the US and Canada should be considered as being covered by UL's Listing and Follow-Up Service meeting the appropriate requirements for US and Canada.**

The UL Listing Mark for the US and Canada generally includes: the UL in a circle symbol with "C" and "US" identifiers:  the word "LISTED"; a control number (may be alphanumeric) assigned by UL; and the product category name (product identifier) as indicated in the appropriate UL Directory.

**Look for the UL Listing Mark on the product**

Issued by: *Walter Hofmair*  
**Walter Hofmair, Senior Project Engineer**

UL International Germany GmbH

Any information and documentation provided to you involving UL Mark services are provided on behalf of Underwriters Laboratories Inc.  
For questions in Germany, you may call +49 0 6102 369 0.

Reviewed by: *Manfred Müller*  
**Manfred Müller, Senior Project Engineer**

UL International Germany GmbH

# Certificate of Compliance

Certificate Number 090905 - E215059  
Report Reference E215059, June 17th, 2005  
Issue Date 2005 September 9

Page 2 of 2



This is to verify that representative samples of the product as specified on this certificate were tested according to the current UL, cUL requirements.

## Standards:

UL 61010A-1 - Electrical Equipment for Laboratory Use; Part 1: General Requirements;  
UL 61010A-2-020 -Electrical Equipment for Laboratory Use; Part 2: Particular Requirements for Laboratory Centrifuges;  
CSA C22.2 No. 1010.1-92 - Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use; Part 1: General Requirements;  
**CSA C22.2 No. 1010.2.020-92, CSA C22.2 No. 1010.2.020A-97 - Part 2: Particular Requirements for Laboratory Centrifuges**  
**UL 61010-1 Electrical Equipment for Measurement, Control, and Laboratory Use; Part 1: General Requirements**  
**CSA22.2 No. 61010-1 - Electrical Equipment for Measurement, Control, and Laboratory Equipment - Part 1: General Requirements**

Issued by:

**Walter Hofmair, Senior Project Engineer**

UL International Germany GmbH

Any information and documentation provided to you involving UL Mark services are provided on behalf of Underwriters Laboratories Inc.UL International Germany GmbH

Reviewed by:

**Manfred Müller, Senior Project Engineer**

UL International Germany GmbH

Any information and documentation provided to you involving UL Mark services are provided on behalf of Underwriters Laboratories Inc.UL International Germany GmbH