

Key features

- Simple, yet exceptionally powerful and versatile centrifuge
- 15000 RPM 15300 RCF
- Patent pending digital motor driver SkySpin(tm) that insures extraordinary performance, low noise, low heating and smooth operation.
- Hight tech digital imbalance sensor with great tolerance.
- Advanced new chassis insure super-smooth run even when slightly imbalanced.
- Secure lid lock.
- 5 rotor break levels.
- Intuitive and clear interface with large color display
- Display preferences include the following:
 - Language
 - Brightness
 - Sleep mode
 - Sound volume
- Short Spin button

Centrifuge **CM50** 3rd generation



Interface



Custom-designed Fugamix interface offers outstanding clarity and comfort.




Using the sound, brightness and language options everyone can make the machine more personal and adequate the current operation environment.

Rotor



Universal high speed rotor for 12 test tubes 0.2 - 2 ml.
Patent pending cover cap insures the optimum air flow and incredibly low sample heating.



Technical specifications

Centrifugation speed	RPM	from 1000 up to 15000
Maximum centrifugal force	RCF	15300
Timer	min	from 0.1 up to 99
Dimensions (L x W x H)	mm	200x180x145
Power supply adapter		24 VDC, 5A
Weight	Kg	3.1