

Minifuge Smart 6

New Product Update: Convenience in the lab



FEATURES

- One spin action simply by closing the lid
- Safety stop by opening the lid
- Selectable operation modes between short spin and continuous 15 minutes
- Interchangeable rotors for 6 micro tubes and 2 PCR strips
- Electro sensitive non-contacting lid sensor

Personal centrifuge for your life science applications.

Special Price \$275

Technical Specifications

Max. RPM	6,500 RPM
Max. RCF	2,234 x g (1.5/2.0 mL Tube) / 2,182 x g (PCR Strip)
Capacity	6 x 1.5/2.0 mL, 2 x 8 Tube PCR Strip
ACC Time to Max Speed	3 secs
DEC Time from Max Speed	4 secs
Operation Mode	Selectable, Short Spin (10 sec) or Continuous (15 min)
Dimension	(W x D x H) 150 x 150 x 90 mm
Weight	650 g
Power Supply	C 100~240V 50/60Hz 2 A
Cat no.	AC Smart 6



The world's best medical technology — in the lab or at the point of care

Order online at www.pocdsscientific.com.au or call **1800 640 075**