


## Technical Data Sheet

Eppendorf Deepwell Plate 96/2000  $\mu$ L

English (EN)



### Plate dimensions [mm]

Length			Width			Height		Stack height
l	l1	l2	w	w1	w2	h	fh	sh (without lid)
127.8	124.5	124.0	85.5	82.2	81.5	44.1	6.1	86.8

### Well dimensions [mm]

Well dimensions [mm]				Bottom [°]	Positions of well center [mm]				Bottom space [mm]	
Depth		Well width		Angle						
h1	h2	d	d1	a	P1	P2	P3	P4	b	b1
38.1	3.1	8.2	7.5	90	14.4	9.0	11.2	9.0	1.5	1.4

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf AG · 22331 Hamburg · Germany

[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)