

Eppendorf Certificate

Certificate of Quality – Quality Assurance of epT.I.P.S.[®], epT.I.P.S.[®] LoRetention, ep Dualfilter T.I.P.S.[®], ep Dualfilter T.I.P.S.[®] LoRetention and ep Dualfilter T.I.P.S.[®] SealMax

All epT.I.P.S. pipette tips and all ep Dualfilter T.I.P.S. filter tips are single-use plastic tips which, in combination with all Eppendorf manual and electronic air-cushion pipettes, form a perfectly functioning pipetting system.

Quality Assurance

Eppendorf hereby confirms that all data for random error and systematic error according to EN ISO 8655 of Eppendorf pipettes are valid only when original Eppendorf pipette tips are used.

All variants of epT.I.P.S. pipette tips and ep Dualfilter T.I.P.S. filter tips are manufactured in a controlled area and are subject to regular quality inspections within a statistical process control.

The Eppendorf Quality Management System is certified according to the international standards ISO 9001 and ISO 13485 plus 14001.

Composition

epT.I.P.S. pipette tips and ep Dualfilter T.I.P.S. filter tips are made of virgin polypropylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into the materials Eppendorf uses for the production of pipette tips:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DiHEMDA))
- Plasticizers (softeners)

Eppendorf Certificate

The raw materials used for epT.I.P.S. pipette tips and ep Dualfilter T.I.P.S. filter tips comply with FDA regulations 21CFR § 177.1520 "Olefin Polymers" and 21CFR § 178.2010 "Antioxidants and Stabilizers for Polymers" for food containers. Eppendorf does not use slip agents, biocides and plasticizers in the production process of epT.I.P.S. pipette tips and ep Dualfilter T.I.P.S. filter tips.

epT.I.P.S. pipette tips and ep Dualfilter T.I.P.S. filter tips and all product components do not contain:

- Metallic dyes
- Latex
- Glucan
- Cellulose

Hamburg, April 2016

Page 2 of 2



Brigitte Klose
Global Product Manager
Liquid Handling Consumables/Instruments

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com



Dr. Birgit Schreiber
Vice President Quality Management &
Regulatory Affairs

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified