

Eppendorf Certificate

Certificate of Quality

Membrane Filter, 0.2 µm for Easypet® 3
International Order No.: 4430 606.005

The above mentioned product, distributed by Eppendorf AG, is subject to careful quality inspection prior to release for distribution.

Each lot is tested 100 % during manufacture for membrane and housing integrity by a leakage test under automatic conditions and by visual inspection.

Each lot is subject to further tests to assure that this product is only released, if it fulfills the requirements and meets the specifications of the technical data.

Sterility is confirmed according to ISO 11135-1.

Hamburg, May 2014



Brigitte Klose
Global Product Manager Consumables



Dr. Birgit Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified