



Ripette® pro

THE HANDY ONE patent pending

Flexible: With the Ripette® pro you can work flexibly with 12 practical tips sizes, 10 volume adjustments and 120 program steps.

Clear: Volume adjustments and dosing vision panel on the Ripette pro.

Handy: A perfect ergonomic design makes easy and comfortable use possible. With the easy-moving dosing button, in combination with the low weight, a long dosing series can be conveniently dispensed.

Compatible: The Ripette® pro is compatible not only with PD-tips from Brand but also with combitips and combitips plus from Eppendorf and has been precisely aligned to Ritips professional from Ritter.



Independent: In addition to the handy design, the Ripette® pro has been constructed with a continuous-action mechanism. The Ripette pro is maintenance-free and does not need any additional resources such as a power supply or batteries.

| Tip size | Program number | | | | | | | | | |
|----------|----------------|------|------|------|------|------|------|-------|------|-------|
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
| 0,05 ml | 1 | 1,5 | 2 | 2,5 | 3 | 3,5 | 4 | 4,5 | 5 | 5,5 |
| 0,10 ml | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 |
| 0,20 ml | 4 | 6 | 8 | 10 | 12 | 14 | 16 | 18 | 20 | 22 |
| 0,50 ml | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 |
| 1,00 ml | 20 | 30 | 40 | 50 | 60 | 70 | 80 | 90 | 100 | 110 |
| 1,25 ml | 25 | 37,5 | 50 | 62,5 | 75 | 87,5 | 100 | 112,5 | 125 | 137,5 |
| 2,50 ml | 50 | 75 | 100 | 125 | 150 | 175 | 200 | 225 | 250 | 275 |
| 5,00 ml | 100 | 150 | 200 | 250 | 300 | 350 | 400 | 450 | 500 | 550 |
| 10,00 ml | 200 | 300 | 400 | 500 | 600 | 700 | 800 | 900 | 1000 | 1100 |
| 12,50 ml | 250 | 375 | 500 | 625 | 750 | 875 | 1000 | 1125 | 1250 | 1375 |
| 25,00 ml | 500 | 750 | 1000 | 1250 | 1500 | 1750 | 2000 | 2250 | 2500 | 2750 |
| 50,00 ml | 1000 | 1500 | 2000 | 2500 | 3000 | 3500 | 4000 | 4500 | 5000 | 5500 |
| | 48 | 31 | 22 | 18 | 14 | 12 | 10 | 9 | 8 | 7 |

Steps

Contamination-free: Ripette® pro and Ritips®

The Ripette® pro is designed for dispensing contaminated fluids and volatile liquids. Due to the exclusive acceptance of dispenser tips with plunger i.e. Ritips and other compatible tips, the Ripette® pro is secure from contamination by ascending aerosols in the tip and guarantees long-term secure use even with samples of high fluidity.

- Due to its low weight the Ripette® pro is most-suitable for long dispensing series.
- Ergonomic design which allows a one-hand-operation of volume selection, loading and dispensing of liquids.
- The dosing table shows all dosing possibilities at a glance for choosing the most-suitable tip for your application.
- The Ripette® pro is reliable and maintenance-free due to its robust design and a minimum of wear parts.
- By using selected premium materials the Ripette® pro shows a high resistance towards a broad range of chemicals.
- Every device will be tested in terms of functionality and precision after production. Only after the requirements of DIN EN ISO 8655 have been fulfilled, a serial number will be issued.

| PARTS No. | Description | Qty | Price (ex. GST) |
|--------------|--------------|-----|-----------------|
| 42020 - 0000 | Ripette® Pro | 1 | \$ 397.00 |

The world's best medical technology — in the lab or at the point of care

Order online at www.pocdscientific.com.au or call **1800 640 075**

Ritips® Standard

The classical range of dispenser tips



| PART No. | Description | Qty | Price (ex. GST) |
|--------------|---------------------------|-----|-----------------|
| 40002 - 0000 | Ritips® 0.05 mL | 100 | \$127 |
| 40001 - 0001 | Ritips® 0.5 mL | 100 | \$90 |
| 40002 - 0001 | Ritips® 1.25 mL | 100 | \$90 |
| 40001 - 0003 | Ritips® 2.5 mL | 100 | \$90 |
| 40001 - 0004 | Ritips® 5.0 mL | 100 | \$90 |
| 40002 - 0002 | Ritips® 12.5 mL | 100 | \$90 |
| 40001 - 0008 | Ritips® 25.0 mL | 100 | \$109 |
| 40001 - 0009 | Ritips® 50.0 mL | 100 | \$109 |
| 40002 - 0010 | Ritips® 0.05 mL bioclean® | 100 | \$233 |
| 40001 - 0010 | Ritips® 0.5 mL bioclean® | 100 | \$161 |
| 40002 - 0011 | Ritips® 1.25 mL bioclean® | 100 | \$161 |
| 40001 - 0012 | Ritips® 2.5 mL bioclean® | 100 | \$161 |
| 40001 - 0013 | Ritips® 5.0 mL bioclean® | 100 | \$161 |
| 40002 - 0012 | Ritips® 12.5 mL bioclean® | 100 | \$161 |
| 40001 - 0015 | Ritips® 25.0 mL bioclean® | 100 | \$194 |
| 40001 - 0016 | Ritips® 50.0 mL bioclean® | 100 | \$194 |

| | | | |
|--------------|-----------------------------|---|-----|
| 43001 - 0847 | Ritips® adapter 25mL, red | 1 | \$7 |
| 43001 - 0848 | Ritips® adapter 50 mL, grey | 1 | \$7 |

Precision dispenser tips standard and bioclean®, suitable for: Ritter Ripette, Ritter Ripette genX®, Multipette® 4780, Minilab 100/101, STEPMATE, EasyStep, Distriman®, HandyStep®, Multipette® plus1, HandyStep® electronic.

1) when using the Multipette® plus the screen will not display.

The following trademarks belong to; Multipette® to Eppendorf AG, HandyStep® to Brand GmbH + Co. KG, Distriman® to Gilson Inc., Bioclean® to Ritter GmbH

INTRODUCTORY SPECIAL

Purchase 1
Ripette® Pro
repetitive
dispenser

and get

2 cases of
RITIPS® STANDARD
of your choice

Please contact your
POCD Scientific
Sales Representative for details!

Call 1800 640 075

The world's best medical technology — in the lab or at the point of care

Order online at www.pocdscientific.com.au or call 1800 640 075

POCD Scientific | Unit B, 19-21 Loyalty Road, North Rocks NSW 2151 | info@pocdscientific.com.au